# THE ALPHABOND REVIEW

#### REVOLUTIONARY PRODUCTS FOR YOUR DENTAL NEEDS

www.alphabond.com.au ISSUE 08 2016



#### GIVE A HIGH LUSTRE SHINE TO ALL CERAMIC RESTORATIONS

Zircon-Brite is a solid polish paste that gives all types of crown and bridge restorative materials a high lustre in a matter of seconds. It is the ultimate polish for Lithium Disilicate, Zirconia, Layered Porcelains, Lucite-based ceramic, CAD/CAM restorations and Composites. Eliminate the need to re-glaze a restoration following adjustments, thereby saving time by having to progress through a re-glazing cycle in a furnace, and exposing the restoration to yet another high-temperature firing cycle. (Order code: HS-DV-ZIRCON20)

#### YOU'LL BE SAYING WOW WITH WOW! SHINE BRIGHT LIKE A



Specifically designed for zirconia, DVA's WOW polishing paste seals and polishes your zirconia restorations to a high degree lustre. In a matter of seconds you'll see how WOW brings out a shine without eradicating special effects on your restoration. (Order code: HS-DV-WOW20)

#### HELP WHEN YOU NEED AN EXTRA HAND

Use TacGel to hold veneers, crowns and small castings in place during fabrication and finishing. With the small dispensing tip, you can use minute amounts with control and precision. Consisting of a specially formulated, water-soluble adhesive, the gel will immediately stabilize the restoration and mating surface. Drying in a short time, it produces a light adhesive bond of the restoration to the surface and then all you need do is evaluate the fit, form or shade of the restoration and make the necessary changes.

Once you're finished, simply remove the TacGel, and

Once you're finished, simply remove the TacGel, and any residual material will burnout leaving no trace of contamination. (Order code: HS-DV-TAGEL)



#### SHINE BRIGHT LIKE A DIAMOND WITH KENDA

Ideal for high shine polishing of flexible partials as well as full dentures. This universal high-lustre polishing paste is also perfect for composites, gold and chrome cobalt. In a stick applicator with approx. 75 ml contents. (Order code: KENDA)



#### A 'TAK'-TIC YOU'LL APPRECIATE




Clinically preferred, Tak's hydroplastic is a custom formulated plastic that softens in hot water. Comfortably handled, staying in place for extended working times.

Hydroplastic hardens to a tough, unbreakable plastic following a five minute bench set, or instantly in cold water. There is no mixing or polymerizing and therefore no odour and no chemical hazard. Capturing minute detail with extraordinary accuracy, Hydroplastic is ideal for customised impressions and bleaching trays. Think of it as wax, where over extensions may be trimmed with scissors, and adjustments made with a flame. Tak's pre-rolled Hydroplastic is excellent for border molding, and available in packets of nine sticks.

(Order code: TAK-STICK9)







### PRIM & PROPER PRIMER

A specialised primer for all aesthetic indirect restorations. GC's Ceramic Primer II is designed to produced a high and durable bond strength to aesthetic materials. No need for refrigeration, you can store the bottle at room temperature. Wave goodbye to mixing, thanks to a single bottle delivery, which allows for fast and easy application. Simply apply and dry.

(Order code: GC-008551)



## RELINE WITH A SOFT TOUCH

Flexor CC is cold-curing and fulfills all demands to make it much more comfortable to wear. It is an automixing, permanently soft direct relining material, based on addition-curing silicone and supplied in cartridges (Automix). The set comes complete with a 50ml silicone cartridge, primer, lacquer catalyst, lacquer basis, polisher, brush, mixing tips and mixing pad. Refill cartridges also available. (Order code: SC250030)



With zirconia's growing popularity, the hardness and sensitivity to heat has been an

issue, but with MPF's new ZMax range you'll cut and abrade with ease. Containing a very high content of fine diamond particles, coupled with special elements and binders, they provide an exceptional lifespan of 2-3 times more than any product on the market today. ZMax requires very little pressure from the handpiece, thus generating less heat to the zirconia substructure.

(Order codes: HP120-0001 and HP120-0001)



#### STERILE EXAMINATION KITS

These single-use instruments are designed to assist in reducing both crosscontamination risks reprocessing costs.

Kits feature high quality, precision-made dental instruments with stainless steel tips and ergonomic latex free handles. All kits are CE and TGA certified and arrive at your surgery or laboratory individually packaged, sterile and ready for single-use with a 2.5 year expiry. Kits available in 10 pack sets or 100 pack sets. (Order code: W9444151)



# PRECIOUS

Repair minor miscasts on precious ceramic alloys, close surface pits, porosity, and small holes prior to porcelain application. Processed in the porcelain furnace, Nu-Met results in a finished metal with bonding characteristics similar to that of a ceramic alloy. Available for Platinium or Gold.

(Order codes: NUMETP or NUMETG) )

dental

