

THE ALPHABOND REVIEW

REVOLUTIONARY PRODUCTS FOR YOUR DENTAL NEEDS

www.alphabond.com.au

ISSUE 12

2017

PROSTHESIS@BREDENT

INTRAORAL & LABORATORY DENTURE REPAIR

Repair dentures quickly with a *Qu-resin*! A quick-setting, self-curing denture repair resin available in two colours (pink and dentine). Based on diacrylate for direct and indirect use. Each assortment set contains: 1 x *Qu-resin* 50ml cartridge, 1 x 10ml Qu-connector, 12 x blue size 1 mixing cannulas.

REF: BR54001161 (*Pink*) or BR54001162 (*Dentine*)



GO AHEAD.. COAT THAT CrCo!

Ropak Compact Opaquer UV is an alternative for conveniently coating Cobalt Chrome frameworks. Simply apply with the integrated brush directly onto the clean metal surface to provide the future acrylic area with an enhanced visual appeal. Available in a 20ml bottle and is ready-to-use.

REF: BR54000133



TOOTH COLOURED OPAQUER

Compact Opaquer tooth-coloured UV enhances esthetics in the area of acrylic teeth. Especially suited for the anterior area, it features a fine coat which is easily applied to the framework, to enhance your restorations visual appeal. Available in a 10ml bottle and is ready-to-use.

REF: BR54000105



SEPARATE & SHINE

Isoplast is a separating liquid for plaster-resin which is suitable for all cold-curing and heat-curing resins.

Isoplast is a separating liquid for plaster-resin which is suitable for all cold-curing and heat-curing resins. Using Isoplast will allow you to obtain a brilliant mirror-finish on resin surfaces. Available in a 750 ml bottle.

REF: BR54001019



DIACRYL ACRYLIC GRINDING INSTRUMENTS



A diamond grinding tool with a recommended speed of 10-15,000 rpm. The fine, pointed flame design allows filigree finishing of alveolar and papillary bases.

REF: BR34001050



A diamond grinding tool with an extra coarse diamond grit and large diameter, it grinds aggressively to create an optimum ground surface. Recommended speed: 10,000 rpm.

REF: BR34001030



A diamond grinding tool for use in lingual and palatal areas, specifically for grinding large papillae and root attachments or lingual bars. Recommended speed of 10-15,000 rpm.

REF: BR34001040



This abrasive rubber instrument creates smooth surfaces on acrylic dentures instead of using sandpaper. With its smooth surface, the denture can be prepared for optimum polishing in minimum time. Recommended speed of 8-10,000 rpm.

REF: BR34000900



ALPHABOND dental

www.alphabond.com.au Ph: (02) 9417 6660 E-mail: info@alphabond.com.au



Like us on Facebook

CROWN & BRIDGE



INTRODUCING BIONICS HIGH-PERFORMANCE POLYMER (BioHPP)



BioHPP is a high-performance polymer (PEEK) material with ceramic fillers. PEEK has been used as implant material in human medicine for over 30 years. Cheaper than gold, BioHPP is easier and better to process in dental laboratories than non-precious metals, titanium or ceramic. Made of extraordinarily elastic material, close to that of bone, it's capable of compensating for the torsion of the bone, particularly in larger implant work. Available as a breCAM. circular blank for processing in CAD/CAM milling machines such as the ZirkonZahn system and others. Available as special order-in items, these blocks are available in diameters of 98.5 mm and 95 mm with various heights and shades. These can be milled with specialised breCAM. cutters, engineered to avoid overheating and developed for the dry processing of BioHPP (PMMA/PEEK/wax) and CAD/CAM composites.

DENTAL ATTACHMENT ADHESIVE

DTK-Adhesive is dual-hardening composite adhesive specially designed for the fixation of dental attachment elements.

Contains 1 x 5g Base paste, 1 x 5g Catalyst paste, 1 x Mixing block, 1 x Spatula and 10 x Application cannulas (black). REF: BR54001185



SUPER-HARD DIE STONE

If you're looking for a synthetic super-hard class IV stone with distinctive thixotropy and improved flow, Exakto-Rock S is exactly what you'll want! - Free from formaldehyde ensures your safety when processing and can be used for holistic dentures with ease. With Exakto-Rock S, reduced expansion is completed after 2 hours and is only 0.08%, thus creating highly accurate impressions with precision-fit restorations. Special dye pigments optimize Exakto-Rock S' reflection of light, making it perfectly suited for scanning in your CAD/CAM system. Available in brown and ivory in quantities of 2kg, 10kg and 20kg. REF: BR5700S_

AIRAQUA TURBINE HANDPIECE

A handy, compact unit with a light-weight handpiece for processing of high performance ceramics, press & metal ceramics. Featuring an air/water spray to avoid overheating, and a smaller rotor for comfort and better vision. REF: BR11001460



PRIMED FOR ADHESION

Achieving adhesion between ceramics and composites can be tricky, using K-Primer ensures adhesion of veneering ceramic, press ceramic, lithium di-silicate, and CAD blanks to your composites. Not for intra-oral use, but it can also be used for repairing ceramic veneers. REF: BRAPK25003



ALPHABOND dental

www.alphabond.com.au Ph: (02) 9417 6660 E-mail: info@alphabond.com.au



Like us on
Facebook